



**INFINITT**  
North America

# G3 End User Guide

---





**ANNA JQUES HOSPITAL**

## Contents

### Section:

- 1- Installing InfinitT G3 PACS on a local PC
- 2- Basic Viewer Functions
- 3- A closer view of tools and buttons
- 4- Removing the G3 viewer from a PC
- 5- Tech tips
  - A- Vista
  - B- Internet explorer setting

## Note

- For windows operating systems only 
- Use Windows Internet Explorer only 
- G3 Viewer Will not operate using FireFox

## Minimum PC requirements

- Operating System- Windows XP or greater
- 1 Gig of RAM Minimum
- 40 gig of free hard drive space
- Monitor resolution- 1024x768 minimum
- Internet Explorer- Version 7 or greater

**Having problems with the PACS?**

**Please call the:**

**Anna Jaques PACS Hotline**

**978-463-1353**

# 1-Installing Infinitt G3 PACS on a Local PC

Open Internet Explorer:

Type in the URL (Address):

<https://starpacs.ajh.org>


Or “Click” on your existing Desktop Icon




**\*\*Assure that the PC is logged in as a Windows Administrator**

Usually your PC's Pop Up Blocker will block the Installation.


Disable the Pop Up Blocker as seen below for the above site ONLY:



Left Click on the Yellow Pop Up blocker bar and choose “Always Allow”.  
(This will only need to be done once)




**Would you like to allow . . . ? Click “YES”**



**Enter User ID and Password**


**Click “Login”**

**Next you will likely notice an Active X security warning:**




**Click “OK”**


**Do you want to continue running scripts on this page? Click “YES”**



**As the system attempts to load the PACS Viewer you will see:**



**Click on the yellow bar. Then choose “Install ActiveX Control” from menu**



**Finally, click “Install”**



**If asked whether to add this URL to “Trusted Sites,” choose Yes**

**Note: You may need to log in one more time by typing in the URL into Internet Explorer before the Infinitt Icon appears on your desktop!**



**After it becomes available, you can simply double click the Icon to Login**


## 2-Basic Viewer Functions

1. Click on the Infinitt PACS logo (to launch the login screen)

2. Enter in your User Name / Password

3. Make sure that **Anna Jaques Hospital** is selected in the location field in the second row at the top. If you see **Amesbury** make sure you use the drop down arrow ▼ and select Anna Jaques Hospital before proceeding.

4. Under "Patient Name", type in the last name only of the patient you are looking for. (You can also search by Patient ID, or any other information you know about the exam.) **Note: Do not use commas**



5. Double click on the exam to load the images.

6. When the exam loads you have tools you can use:

(Hovering over a tool will reveal its function)

Toggle between Open Studies

Go back to the worklist without closing

See the Report

Close the exam(s)

Zoom


Window / Level

Measuring Tools (Ruler, ROI, Pixel Lens, Spine Label, Scoliosis)

Reset Changes

Minimize PACS

Exit PACS




Double Click over an image to maximize it to full screen. Use the wheel on your mouse to scroll through the images.



# Compare Studies

With the exam loaded, change the virtual monitor layout to 1x2 (1 row x 2 columns).



In the virtual monitor toolbar where you'd like to load another study, click the related exams drop down box (pictured) and choose the exam you'd like to view. (The default choice in the drop down is Exam Original, meaning the original exam you loaded)



**notes:**

- 1.) When you click on an image for a study in the VIEWER, the thumbnails will display for that study in the THUMBNAIL bar, allowing you to drag a different sequence onto the VIEWER.)
- 2.) You can use a different layout (like 2x2) to load 4 different studies at the same time using the above steps. However, the more studies at a time you load, the smaller the images will be.)

## ER note

**ER & RAD: (Red)- Both left findings with a Disagreement**

**ER & RAD: (Blue)- Both left findings and agreed**


**RAD: (Yellow)- Only Radiology left findings (Wet Read)**

**ER & RAD: (Purple) Both left findings, but the radiologists did not click agree or disagree**

**ER: (Green)- Only the ER has left findings**




### 3-A closer view of tools and buttons in the Viewer



## 4-Removing the G3 viewer from a PC


### Step-1

Delete viewer folder “INFINITT\_AJH” from hard drive. (My computer> c: )



### Step-2

“Remove” the “Smart Update control” file. (C: > Windows> Downloaded program files)




## 5-Tech Tips


A- Not able to load using Vista:

Most likely the User account Control will need to be turned off. Follow Below

➤ Go to Control panel > User Accounts. "Click" Turn User Account Control on or off



➤ Uncheck, Use User Account Control. Then "OK"




## B- Internet Explorer setting


- On a rare occasion the user logged in is NOT the user in the Viewer.



From your desktop Right Click on Internet explorer and select **“Properties”**



In the middle of the Properties screen **“Click”** on **“Settings”**



select **“Every visit to the page”** Then **“Click”** OK